

When it comes to prescription misuse...

Put a Lid on It!


Prescription Misuse Prevention Kits Product Information

What Gets Measured, Improves!


It's a Stopwatch!



Perfect for
Seniors!



No Alarms!
Nothing to set!



Wide Mouth Kit
fits bottles inside.

Automatically starts
when cap is closed!



Why Choose Prescription Misuse Prevention Kits?

TimerCap's Prescription Misuse Prevention Kits are the perfect tool to assist in preventing medication abuse and mistakes. Its patented stopwatch design makes it the most intuitive and simplest way to protect your families and home. Track the time a pill bottle was last opened, down to the minute! It is as easy as closing the cap to start the timer and opening the cap to stop the timer. There are **no alarms or programs to set, yes it is that simple!** Families and communities can now have peace of mind in knowing how easy it can be to stay on top of taking your medications as prescribe, to detect unexpected openings and to deter others from experimenting.

TimerCap
Prescription Misuse Prevention Kit
Automatically shows the last time opened to avoid medication mishaps.

TimerCap
(800) 557-4072
TAKE ONLY AS PRESCRIBED™

TimerCap
Prescription Misuse Prevention Tips

Patient Safety

- Measure** Dosage and time of medication to avoid accidental overdose.
- Monitor** Your abilities before driving to avoid DUI mistakes.
- Manage** Medication in containers that have label and safety information.

Home Safety

- Detect** Unwanted openings by seeing the exact time it was last opened.
- Deter** Misuse and theft of medication.
- Dispose** Safety of unneeded medication.

Making TimerCap Child Resistant (Ships EZ Twist)


- Put ring on top of the TimerCap.
- With CR Ring on, turn upside down on a hard surface.
- Click. Click.
- Press firmly on bottle bottom until you hear a double click.
- When CR Ring is properly CR Ring will not pull off.

Changing the TimerCap Battery (Battery Lasts 12-18 Months)

- Remove bottom disk and Spring.
- Remove a screw holding battery.
- Replace LR-41 battery slide.
- Snap bottom disk with Spring back in place and attach screw. Done!

Date: _____ Medication: _____

TimerCap
TimerCap.com
(800) 557-4072



TimerCap Mission

One in four teens admit to abusing or experimenting with prescription drugs; of those, 42% of teens get their drugs from their own families home. Many parents mistakenly think that prescription drugs are safer than street drugs. This idea has led to about 15,000 people fatally overdosing on prescription painkillers every year. That is more than cocaine and heroin overdose deaths combined.

At an early age Larry Tersky, TimerCap's CEO, experienced first hand how prescription drugs can affect a family. Growing up witnessing his own mother's addiction to prescription pain killers until her untimely passing, combined with shocking abuse statistics, is what drove Larry to find TimerCap. TimerCap now takes pride in saying we continue to collaborate everyday with anti-drug coalitions, healthcare groups, counties, state, parents, police, schools, youth, and media to provide tools for healthier drug free communities.

Larry Tersky enjoys to lecture in a fun and very informative way on a variety of topics such as communication skills, marketing, prescription misuse prevention, and more. Larry has proven and is determined to continue to associate TimerCap as a first line tool of defense against misuse of medication and to make a difference in stopping prescription abuse. No matter if you take medication for chronic issues such as: pain pills, depressants, opioids, Alzheimer, diabetes, ADHD, HIV, PrEP, or recovery from a procedure, be part of the solution with TimerCap and help change the conversation. Let us join together to bring safety and peace of mind to families and communities.

Standard Style Prescription Misuse Prevention Kits

TimerCap Benefits and Product Features



Color coded for Easy Identification
Easily identify your medication by color for your convenience.




Fits a Large Range
With the very popular 38-400 threaded bottle, these easy and simple caps will fit most large, pharmacy and vitamin bottles.



It's a Stopwatch
It is not a clock, so no alarms, no programs to set! The patented pressure switch automatically starts when the stopwatch cap is closed and stops when cap is opened.



Know, Don't Guess!
Never forget when you last took your medication. The automatic stopwatch feature counts up for you.



What Gets Measured, Improves!



Screen shows time elapsed. NOT TIME OF DAY!



Screen shows: SEC 1, SEC 2 then 10, 11, for the first minute.



Continues to count minutes and hours since last opened, until 99:59, then it will go blank.



Safety First
TimerCap ships EZ-Twist, and it's great for seniors. It also contains a caution warning on each side in both English and Spanish.



Child Resistant Ring
TimerCap comes with a child resistant ring (CRC) that snaps on to protect your children.



Detect and Deter
Due to 52% of misuse starting at home, it is vital to have a tool that can detect, deter, and prevent theft or openings by seeing the time it was last opened.

Wide Mouth Style Prescription Misuse Prevention Kits TimerCap Benefits and Product Features



Color coded for Easy Identification

Easily identify your medication by color for your convenience.



Fits Standard Rx Bottles

Wide Mouth's fit original prescription bottles that are within its width and height limits of 1^{3/4}" wide and 3^{1/4}" tall.




It's a Stopwatch

It is not a clock, so no alarms, no programs to set! The patented pressure switch automatically starts when the stopwatch cap is closed and stops when cap is opened.



Don't Guess, Know!

Never forget when you last took your medication. The automatic stopwatch feature counts up for you.



What Gets Measured, Improves!

HOURS:MINUTES

02:00

Time Since Last Opened

Screen shows time elapsed. NOT TIME OF DAY!

HOURS:MINUTES

SEC 1

Time Since Last Opened

Screen shows: SEC 1, SEC 2 then 10, 11, for the first minute.

HOURS:MINUTES

99:59

Time Since Last Opened

Continues to count minutes and hours since last opened, until 99:59, then it will go blank.



Safety First

Wide Mouth Caps are available in **EZ Twist (Red Cap)** and in a **Child Resistant Cap (White Cap)** for your child's safety.



Warnings for EZ-Twist off

TimerCap's EZ-Twist has a caution warning on each side, in both English and Spanish, making it clear that child resistant (CRC) ring is **NOT** on.



Detect and Deter

Due to 52% of misuse starting at home, it is vital to have a tool that can detect, deter, and prevent theft or openings by seeing the time it was last opened.

Standard and Wide Mouth Specifications Sizes and What is Included In Your Kit




Standard Style Kit
(Convertible Cap with child resistant ring)



Wide Mouth Style Kit
(Choose between EZ Twist or CRC-Child Resistant)

Cover Sheet - Front




The cover sheet will give you images of what is included in the kit of your choice. It also will display e-mail, phone, or web contact information.

Medication Tracker




The Medication Tracker form will allow you to keep track of when your pills were taken, how many were taken/remaining, and how one felt.

Prevention Tips



Prevention tips are safety tips or guided reminders of how prescribed medication should be monitored in ones home in order to stay safe.

Instructions - Back



The back instructions will give information on how certain features of specific kits work and how the stopwatch feature counts time.

Side Label



The side label is provided in case one does not want, or can not, peel off the original prescribed label. Just fill in the date it was filled and the medication.

TimerCap Customizable Prescription Misuse Prevention Kits

Design Areas and What is Included in Your Custom Kit

Get any colored "D" style or Wide Mouth TimerCap with your choice of design on the label. We also updated our instructions and tracker forms to give our customers more room for their choice of information. Make TimerCaps yours by adding your company logo, contact information, or even a message. Let TimerCap be a part of your change to the community.



Standard Style Kit
(Convertible Cap with child resistant ring)



Wide Mouth Style Kit
(Choose between EZ Twist or CRC-Child Resistant)

Your Company Design Area

[Your_Website.com](#)

Prescription Misuse Prevention Kit

Automatically know the last time opened to avoid medication mishaps.



TimerCap



CR Ring



Bottle



Tips and Tracker

Any questions, complaints, compliments, custom orders, insights, issues or suggestions, please contact us via e-mail at Info@TimerCap.com or visit TimerCap.com. Thank you!



Your Company Design Area

(555) 555-555

**** TAKE ONLY AS PRESCRIBED ****

Patient _____ Medication _____ Dosage _____ Pill Count _____

Date	Time on Cap	Feel (0-10)	(Pills) Taken	(Pills) Remain	Date	Time on Cap	Feel (0-10)	(Pills) Taken	(Pills) Remain

Your Company Design Area

Prescription Misuse Prevention Tips

Patient Safety



Measure
Dosage and time of medication to avoid accidental overdose.



Monitor
Your abilities before driving to avoid DUI fatalities.



Manage
Medication in containers that have label and safety information.

Home Safety



Detect
Unwanted openings by seeing the exact time it was last opened.



Deter
Misuse and theft of medication.



Dispose
Safety of unneeded medication.

Date: _____ Medication: _____

Your Company Design Area

Visit TimerCap.com, call us at (800) 557-4072 or email us at info@timercap.com if you wish to place a customized order.

Let us Join Together to Change the Conversation How to Bring Awareness to Your Community


"I don't have enough staff to put a body everywhere, but I can certainly send TimerCap sets. Everyone loves these!"

Lisa Perry, Ph.D. Director of Prevention Services
Crowley's Ridge Development Council, Inc.

You are Not Alone, TimerCap is Here to Help

Motivating the community to be vigilant about their medication is something coalitions struggle with all the time. TimerCap has dedicated countless hours into working with coalitions to design a turn-key kit that covers the most important issues about preventing prescription misuse. With patient safety in mind, we have created a cost effective tool that drastically increases the number of homes protected. These prescription misuse prevention kits along with our safety tips and medication tacker form, result in some of the best practices that you can use to improve medication vigilance in your community. Typically when coalitions discuss prescription misuse, it starts with statistics leading to limited success.

When starting with statistics like the image below, most parents will give you the glazing eyes along with **"Well, not in my home."** Instead, of starting with statistics of people they don't know or can not relate to, try starting the conversation centered on the person taking the medication and their safety, like the quote below.



"Has anyone in your home ever forgotten if they did or did not take their medication?"

This sparks people to raise their hands in large meetings or it starts a very productive one on one conversation as the answer is almost always **"Yes,"** followed up with a story of who and how it affected the person. Visit TimerCap.com to [download this full script for free](#), and use it at your next big event. Coalitions have seen amazing results using this script and the misuse prevention kits to create awareness and change in their communities.

Visit TimerCap.com to order your 100 pack kit today! Choose between or Standard Convertible Cap, Wide Mouth EZ Twist or our Wide Mouth Child Resistant Prescription Misuse Prevention Kit.

Standard EZ/CRC Cap Kit



Wide Mouth EZ Twist Kit



Wide Mouth (CRC) Kit

